

CS Accessory Case

Equipment for Calibration



Applications

- ✓ Case with all important standard accessories for the calibration of accelerometers
- ✓ Effective calibration of accelerometers in the laboratory
- ✓ Applies for nearly all sine and shock calibration systems (we offer a separate accessories case for CS18 HS)
- ✓ Covers most of the well-known manufacturers (like Endevco, Brüel & Kjaer, Metra, Kistler, PCB, Dytran)

Features

- ✓ Useful items for handling of DUT's before, during and after a calibration
- ✓ Broad variety of mounting adapters for DUT mounting either on Back-To-Back-Sensors or directly on shaker heads
- ✓ Selection of electrical adapters for sensor cables



Components

Quantity	Item
Accessories for the mechanical mounting of sensors	
2	Set of Allen Key (1x imperial; 1x metric)
1	Special adapter set for PR-sensors e.g. Endevco 7264
1	Torque wrench
7	Nuts for torque wrench (9/32"; 1/2"; 3/4"; 1-3/16"; 16mm; 17mm; 22mm)
1	Adjustable wrench
1	Adapter plate with a variety of mounting adapters and test threads
2	Slot screwdrivers
Accessories for the electrical connection of sensors	
12	Adapters for BNC, TNC and Microdot
3	Sensor cables 3m (1x BNC - Microdot, 1x TNC - Microdot, 1x Microdot - Microdot)
Additional additives	
1	Knife
1	Glass bottle for Acetone
1	Duct tape
1	Loctite 401 Super-Glue
1	Gluing wax
1	Small box with vaseline
1	Insulation tape
1	Sandpaper
1	Folding rule
2	Pencils (1x permanent, 1x non-permanent)
10	Self-adhesive cable clips